

## QUICK MEDIC POUCH MK2

Part Number

**Q.M.P.2. 0830/1**

Color

**MULTICAM**

Composition

**100% POLIAMMIDICA**

Other materials

**95% POLIAMMIDICA**

**5% ELASTAN**



UNIQUE SIZE

### WASH INSTRUCTIONS

Wash



Bleaching



Ironing



Dry wash



Drying



Product by: **Brooklin S.r.l.**

Via XI Febbraio,27 - 50053 Empoli (FI) - Tel. 0571 73041/42 - Fax: 0571 72175

Pec: brooklyn@pec.it - e-mail: brooklyn@brooklyn.it

C.C.I.A.A. Firenze 289096 - Trib. 29196 - Pos. M. FI003140 - Part. IVA e Cod. Fisc. it01703650489

Nato Cage AG 941 - Iscritta al n. 609 del Registro Nazionale delle Imprese del Ministero Della Difesa - Presente sul MEPA come Brooklin S.r.l.

## QUICK MEDIC POUCH MK2

### MODEL DESCRIPTIONS

Quick personal medic pouch. Can be taken from the left or the right side. Thi system is made by an outer holding interface and an inner pouch that keep the personal first aid kit. The inner pouch has enough place to keep a complete Medical kit. Tegris interface compatible with Molle or belt.

The structure is made in Cordura 500D, Tegris and an abrasion resistant Schoeller stretch fabric. Totally laser cut. Weight of the pouch 100gr ca.



Part Number  
Q.M.P.2. 0830/1

Product by: **Brooklin S.r.l.**

Via XI Febbraio,27 - 50053 Empoli (FI) - Tel. 0571 73041/42 - Fax: 0571 72175

Pec: brooklyn@pec.it - e-mail: brooklyn@brooklyn.it

C.C.I.A.A. Firenze 289096 - Trib. 29196 - Pos. M. FI003140 - Part. IVA e Cod. Fisc. it01703650489

Nato Cage AG 941 - Iscritta al n. 609 del Registro Nazionale delle Imprese del Ministero Della Difesa - Presente sul MEPA come Brooklin S.r.l.